

## ***Advanced Temperature-Programmed Oxidation (TPO) of Coked Catalysts Using Methanation and FID Detection***

Today, advanced techniques such as TEM, Laser Raman Spectroscopy, EELS,  $^{13}\text{C}$  NMR, and temperature-programmed oxidation (TPO) are widely used to study coked catalysts. Among these, TPO has become one of the most commonly applied methods due to its simplicity and effectiveness.

This *Altamira Note* discusses the use of TPO combined with an innovative detection method developed by Dr. S.C. Fung and Dr. C.A. Querini at Exxon Research and Engineering Company. This approach is straightforward and enables continuous monitoring of the rate of coke oxidation.


### **Experimental**

In this TPO method,  $\text{CO}_2$ —a gas to which flame ionization detectors (FID) are typically insensitive—is converted to  $\text{CH}_4$ , which is easily detected by an FID. The conversion occurs in the presence of a carrier gas containing oxygen.

For these experiments, an **AMI Catalyst Characterization System** equipped with a methanator and FID was used. *Figure 1* shows the system's flow diagram.

The methanator consisted of a small reactor filled with a ruthenium catalyst, positioned downstream of the sample U-tube. When hydrogen passed through the methanator, the Ru catalyst quantitatively hydrogenated  $\text{CO}_2$  to  $\text{CH}_4$ . The FID continuously monitored the rate of  $\text{CH}_4$  formation, providing a real-time measurement of the coke oxidation rate.

A GC column was unnecessary because the FID is insensitive to oxygen and water vapor in the gas stream.



For these experiments, approximately 20 mg of coked catalyst were loaded into the sample cell. A helium carrier gas containing a low concentration of oxygen flowed over the sample at a rate of 20–80 cc/min. The temperature was increased linearly from room temperature until the complete oxidation of all carbon deposits.

The methanator contained approximately 500 mg of **40 wt% Ru/zeolite 13X**. A pure hydrogen stream was injected into the methanator at a flow rate of 22 cc/min. Under these conditions, **CO<sub>2</sub>** was quantitatively converted to **CH<sub>4</sub>**, while the oxygen in the carrier gas was reduced to water. The combined gas stream then flowed directly into the FID.


The FID continuously monitored the methane generation rate, which was equivalent to the rate of coke oxidation.

### **Influence of Oxygen Concentration, Flow Rate, and Methanator Temperature**

Since TPO experiments require an excess of oxygen, it was necessary to evaluate how oxygen concentration affects the hydrogenation of **CO<sub>2</sub>** to **CH<sub>4</sub>** and to establish optimal operating conditions. The effects of flow rate and methanator temperature on the hydrogenation efficiency of the ruthenium catalyst were also investigated, along with the impact of various pretreatments on the Ru catalyst's activity.

To study the influence of oxygen concentration, pulses of **1%, 2%, and 4.26% CO<sub>2</sub>** in helium were introduced into the methanator using a pure helium carrier. Additional CO<sub>2</sub> pulses were introduced using helium carriers containing **0.5%, 1%, and 3% oxygen**.

As shown in *Table 1*, the CO<sub>2</sub> pulses were completely converted to CH<sub>4</sub> except when the oxygen concentration was increased to **3%**. One possible explanation for this behavior is that water formed in the methanator (due to the oxygen present in the carrier gas) reduced the equilibrium conversion of **CO<sub>2</sub> to CH<sub>4</sub>**.



<b>Table 1</b>					
Effect of O <sub>2</sub> Concentration on the Conversion Efficiency of a Methanator—CO <sub>2</sub> Pulse to CH <sub>4</sub>					
% CO <sub>2</sub> in Pulse	% O <sub>2</sub> in He Carrier	% CO <sub>2</sub> Conversion to CH <sub>4</sub>			
		0	0.5	1	3
1		100	100	100	63
2		100	100	100	61
4.26		100	100	100	55
Methanator at 400°C and carrier gas 60 cc/min					

<b>Table 2</b>						
Effect of Oxygen on the Equilibrium Conversion of CO <sub>2</sub> Hydrogenation to Methane						
The influence of O <sub>2</sub> on the CO <sub>2</sub> methanization equilibrium is via water formation:						
$\frac{1}{2} \text{O}_2 + \text{H}_2 \rightleftharpoons \text{H}_2\text{O} \quad (1)$ $\text{CO}_2 + 4 \text{H}_2 \rightleftharpoons \text{CH}_4 + 2 \text{H}_2\text{O} \quad (2)$						
Equilibrium CO <sub>2</sub> Conversion						
%O <sub>2</sub>	1% CO <sub>2</sub> Pulse			4% CO <sub>2</sub> Pulse		
	320°C	400°	430°C	320°C	400°C	430°C
0	100	99.99	99.98	99.99	99.60	98.52
1	100	99.98	99.91	99.99	99.13	96.91
3	99.99	99.85	99.40	99.94	99.26	90.06
10	99.74	86.00	66.45	78.00	49.70	37.32

Higher oxygen concentrations lead to greater water formation, which in turn affects the equilibrium conversion of **CO<sub>2</sub> to CH<sub>4</sub>**. *Table 2* shows that TPO experiments should be conducted with **oxygen concentrations at or below 3%** and at **methanator temperatures below 430°C** to avoid equilibrium limitations. However, the incomplete conversion of CO<sub>2</sub> to CH<sub>4</sub> observed in *Table 1* is *not* attributable to equilibrium constraints.

An alternative explanation is that **water inhibits the methanation activity of the ruthenium catalyst**. To investigate this, experiments were conducted varying three parameters: the oxygen concentration in the carrier gas, the methanator temperature, and the carrier gas flow rate. As shown in *Table 3*, it was necessary to reduce the flow rate when higher oxygen concentrations were used

Additional experiments introduced water directly into the system by saturating the carrier gas at room temperature, producing a **2.6% water concentration in helium**. This water level is comparable to that generated during oxidation in a **1.3% oxygen**

environment. These results confirmed previous findings that **oxygen concentrations should ideally remain below 2%**. For experiments requiring higher oxygen levels, an **oxygen trap** can be installed upstream of the methanator. These traps effectively remove oxygen without affecting the **CO<sub>2</sub> concentration** exiting the sample U-tube.

### Temperature-Programmed Methanation Studies

Temperature-programmed methanation experiments were also performed to evaluate the **stability of the ruthenium catalyst** and to determine the **optimal methanator temperature** for converting CO<sub>2</sub> to CH<sub>4</sub>.

As indicated in *Table 3*, **CO<sub>2</sub> conversion increased** with both rising temperature and decreasing flow rate. This behavior suggests that conversion limitations in the presence of oxygen are primarily **kinetic** in nature. Additionally, the catalyst's activity improved with temperature. However, the methanator temperature should be kept **as low as possible** to minimize the risk of sintering the ruthenium particles, which would permanently reduce catalyst activity.

<b>Table 3</b>			
Oxygen In Carrier Gas Slow the Rate of Methanation via Water Production in the Methanator			
%O <sub>2</sub> in He Carrier	Methanator Temperature°C	%CO <sub>2</sub> Conversion (4.26% CO <sub>2</sub> Pulse)	
		At Carrier Gas Flow Rate	
		60 cc/min	20 cc/min
1	350	82	100
1	400	100	100
3	400	55	100
(2.6% H <sub>2</sub> O)	350	75	100
(2.6% H <sub>2</sub> O)	400	100	

### Effect of Carrier Gas Flow Rate and Catalyst Stability

Experimental results indicated that **FID sensitivity increased linearly** with carrier gas flow rates up to **60 cc/min**. At higher flow rates, the FID response plateaued, suggesting that flow rates above this level do not further improve sensitivity.

An additional important observation was the **deactivation of the ruthenium catalyst** in


the methanator due to sulfur poisoning. This deactivation was caused by sulfur oxides generated during the combustion of sulfur-containing coke deposits. The most effective solution was the installation of a **sulfur oxide trap upstream of the methanator**, which successfully removed sulfur contaminants without affecting the CO<sub>2</sub> concentration.

## Conclusions

These experiments demonstrate that **TPO coupled with methanation and FID detection** is a highly effective technique for monitoring the **carbon oxidation rate of coked catalysts**. By optimizing experimental parameters, complete conversion of **CO<sub>2</sub> or CO to CH<sub>4</sub>** is achievable, even in the presence of oxygen-containing carrier gases.

This method is sensitive enough to detect **carbon concentrations below 0.1%** and can distinguish subtle variations in the coke distribution on catalyst surfaces.

This *Altamira Note* summarizes a presentation delivered by **Dr. S.C. Fung** at an AMI (formally: Altamira's) U.S. User's Meeting. For further details on this TPO methodology, see: S.C. Fung and C.A. Querini, "A Highly Sensitive Detection Method for Temperature-Programmed Oxidation of Coke Deposits: Methanation of CO<sub>2</sub> in the Presence of O<sub>2</sub>," *Journal of Catalysis*, **138**, p. 240 (1992).



Comparison of TPO spectrum obtained with a GC column-TCD and with a methanator-FID. Carrier gas, 1% O<sub>2</sub>/He; carrier flow rate, 60 cc/min; heating rate, 13 K/min.

**Note:** AMI is the only company to integrate a chemisorption analyzer platform with an FID detection system and methanation reactor for advanced TPO studies. This unique configuration enables precise, real-time quantification of coke oxidation rates with unparalleled sensitivity.