

# Overview of Supported Ionic Liquids: Effect of Specific Surface Area and Pore Structure

## 1. Background

Ionic liquids (ILs) are room-temperature molten salts composed entirely of ions—typically an organic cation paired with an inorganic or organic anion. Compared with conventional organic solvents, ILs offer:

- ✓ Broad liquid-phase windows and thermal/physicochemical stability (often stable below ~300 °C),
- ✓ Negligible vapor pressure and nonflammability for cleaner operation, even under high vacuum,
- ✓ Wide solubility ranges that accommodate inorganic/organic compounds and polymers, sometimes with dual roles as medium and catalyst,
- ✓ Favorable electrochemical behavior for electrolytes and reaction media, and;
- ✓ Molecular tunability by varying cation/anion structure.<sup>(1,2)</sup>

Supported ionic liquids (SILs) are functional ILs immobilized on porous carriers—by physical deposition or chemical anchoring—to form thin IL films within a solid matrix. Confining the IL converts a difficult-to-handle liquid into a solid-like composite that combines the chemical selectivity of the IL with the mass-transfer and mechanical advantages of a porous support. This approach mitigates the high viscosity and separation challenges of neat ILs and simplifies recovery and reuse, which is why SILs have become a focus across adsorption and catalysis. SILs show broad utility in gas and liquid adsorption/separation,<sup>(3)</sup> function as catalysts or catalyst supports,<sup>(4)</sup> and provide efficient aqueous-phase removal of heavy-metal ions.<sup>(5)</sup>

For SILs the support's specific surface area (SSA) and pore architecture become critical: high SSA and well-matched pore sizes disperse thin IL films, maximize accessible interfacial sites, and shorten diffusion paths, directly impacting adsorption capacity, selectivity, and kinetics. Routine N<sub>2</sub>/CO<sub>2</sub> physisorption (e.g., with an **AMI Micro 200** or **Sync 400**) to quantify BET SSA and PSD is therefore essential before comparing adsorption performance across varying IL formulations.<sup>(6)</sup>

Common supports feature high specific surface area (SSA) and engineered pore networks—including silica, activated carbon, alumina, MOFs, zeolites, and polymers.<sup>(3,7,8)</sup> Depending on the binding strategy, SILs are described as physically supported (e.g., impregnation) or chemically supported (e.g., covalent bonding, sol-gel routes). In both cases, SSA and pore size distribution (PSD) govern IL film thickness, dispersion, and accessibility. High SSA and well-matched pores deliver uniform IL dispersion, maximize interfacial contact with the fluid phase, and shorten diffusion paths—all of which raise capacity, selectivity, and kinetics. Conversely, excessive IL loading into narrow pores can block pathways and reduce uptake. Quantifying SSA/PSD helps distinguish transport limitations from true site effects and guides optimization of IL loading and support selection.<sup>(6)</sup>

## 2. Representative studies

- ✓ TMGL–SiO<sub>2</sub>: SO<sub>2</sub> capacity comparable to neat IL with improved selectivity; excess loading reduces accessible area and lowers uptake.<sup>(9)</sup>
- ✓ IL–ZIF–IL (core–interlayer–shell): an interlayer ZIF-8 prevents IL aggregation, improving CO<sub>2</sub> capacity and selectivity and maintaining stability over cycling.<sup>(10)</sup>
- ✓ Graphene oxide-based hydrogels: hierarchical porosity with IL functionality affords high heavy-metal capacity and good regeneration.<sup>(11)</sup>

## 3. Case study with AMI Micro 200

Shang et al. prepared a graphene oxide/dicationic ionic-liquid composite (GO-DIL) as a new adsorbent for environmental remediation, with a focus on removing hexavalent chromium, Cr(VI), from water.<sup>(12)</sup> Batch tests showed that pH = 3 yielded the highest Cr(VI) uptake, consistent with maximum electrostatic attraction between protonated sites on GO-DIL and anionic chromium species under acidic conditions. To relate performance to texture, graphene oxide (GO) and GO-DIL were characterized by N<sub>2</sub> adsorption at 77 K using the **AMI Micro 200** static physisorption analyzer (formerly called JW-BK200C). The resulting adsorption–desorption isotherms were used to derive BET specific surface area and pore metrics. These surface and pore data provide the basis for interpreting differences in capacity, kinetics, and regenerability between GO and GO-DIL.

Sample	BET SSA (m <sup>2</sup> /g)	Pore Volume (cm <sup>3</sup> /g)
Graphene oxide (GO)	9.93	0.043
GO-DIL	10.46	0.077


*Table 1: Summary of SSA and pore volume results from Shang et al.<sup>(12)</sup>*

As summarized in Table 1, the BET specific surface area and pore volume of GO-DIL are noticeably higher than those of GO, indicating greater accessible interface and improved mass-transfer pathways. This textural gain aligns with the higher Cr(VI) adsorption capacity observed for GO-DIL and helps explain the performance difference. The materials were further characterized by SEM, XPS, and FTIR, and adsorption–desorption cycling was evaluated. Kinetic fits and Langmuir isotherm analysis point to an electrostatically driven mechanism in which protonated imidazolium sites (N<sup>+</sup>) on GO-DIL attract anionic Cr(VI) species under acidic conditions, enhancing uptake. Overall, the results support GO-DIL as an effective supported-ionic-liquid adsorbent for Cr(VI) and related heavy metals in wastewater and underscore the value of routine N<sub>2</sub>/CO<sub>2</sub> physisorption (BET, pore volume, PSD) for correlating surface/pore metrics with adsorption performance.

## 4. Conclusions

For surface-area and pore-structure measurements that underpin supported-ionic-liquid design, the **AMI Micro 200**, shown in Figure 1, provides the right combination of high-vacuum stability, precise low-pressure control, and robust cryogenic/thermostatic operation (77 K N<sub>2</sub>, 87 K Ar, 195/273 K CO<sub>2</sub>). This enables accurate BET SSA, total pore volume, and PSD across ultralow P/P<sub>0</sub> through near-saturation. Therefore, slight differences in IL loading, film thickness, or support selection can be resolved with confidence.

Automated dose/evacuation routines, tight transducer linearity, and multi-sample throughput improve repeatability for comparative studies, while seamless switching among probe gases links texture (SSA/PSD) directly to capacity, selectivity, and kinetics observed in downstream adsorption tests. In short, the **Micro 200** from AMI delivers the quantitative pore-structure insight needed to optimize SIL formulations for application-specific performance measurements.



- ✓ High precision analysis of specific surface area and pore size distribution
  - ✓ Pore sizes from 0.35 nm – 500 nm
- ✓ Reduced internal dead volume
- ✓ Two integrated degassing ports and two independent analysis ports
- ✓ Optional turbo molecular pump for high vacuum analysis

Figure 1: Highlighted features of **Micro 200** physisorption analyzer from AMI

## 5. References

- (1) Holbrey, J. D. and Seddon, K. R. Ionic liquids. *Clean Products and Processes*, **1999**, *1*, 223-236.
- (2) Kaur, G.; Kumar, H.; Singla, M. Diverse applications of ionic liquids: A comprehensive review. *J. Mol. Liq.* **2022**, *351*, 118556.
- (3) Li, X.; Chen, K.; Guo, R.; Wei, Z. Ionic liquids functionalized MOFs for adsorption. *Chem. Rev.* **2023**, *123*, 10432-10467.
- (4) Wolny, A. and Chrobok, A. Silica-based supported ionic liquid-like phases as heterogeneous catalysts. *Molecules*, **2022**, *27*, 5900.
- (5) Wang, S.; Geng, X.; Zhao, Z.; Zhang, M.; Song, Y.; Sun, K.; Zhang, Q. Ammoniated-driven green synthesis of charged polyoxometalate supported ionic liquids for exceptional heavy metal remediation in actual industrial wastewater. *Water Res.* **2025**, *272*, 122939.
- (6) Rufete-Beneite, M.; Román-Martínez, M. C.; Linares-Solano, A. Insight into the immobilization of ionic liquids on porous carbons. *Carbon*, **2014**, *77*, 947-957.
- (7) Lemus, J.; Palomar, J.; Gilarranz, M. A.; Rodriguez, J. J. Characterization of supported ionic liquid phase (SILP) materials prepared from different supports. *Adsorption*, **2011**, *17*, 561-571.
- (8) Lin, Y.; Wang, F.; Zhang, Z.; Yang, J.; Wei, Y. Polymer-supported ionic liquids: Synthesis, characterization and application in fuel desulfurization. *Fuel*, **2014**, *116*, 273-280.
- (9) Zhang, Z.; Wu, L.; Dong, J.; Li, B.-G.; Zhu, S. Preparation and SO<sub>2</sub> sorption/desorption behavior of an ionic liquid supported on porous silica particles. *Ind. Eng. Chem. Res.* **2009**, *48*, 2142-2148.
- (10) Han, G.; Yu, N.; Liu, D.; Yu, G.; Chen, X.; Zhong, C. Stepped enhancement of CO<sub>2</sub> adsorption and separation in IL-ZIF-IL composites with shell-interlayer-core structure. *AIChE J.* **2020**, *67*, e17112.
- (11) Kong, W.; Chang, M.; Zhang, C.; Liu, X.; He, B.; Ran, J. Preparation of xylan-g-P(AA)-co-AM)/GO nanocomposite hydrogel and its adsorption for heavy metal ions. *Polymers*, **2019**, *11*, 621.
- (12) Shang, J.; Guo, Y.; He, D.; Qu, W.; Tang, Y.; Zhou, L.; Zhu, R. A novel graphene oxide-dicationic ionic liquid composite for Cr(VI) adsorption from aqueous solutions. *J. Hazard. Mater.* **2021**, *416*, 125706.