

# ***Comparison of Static and Dynamic Flow Measurements for Specific Surface Area of Silicon Nitride Powder***

## **1. Background**

With the advancement of science and technology, the demand for novel materials in various industries has increased significantly. Due to the inherent limitations of metallic materials, structural ceramics are gradually replacing them in certain fields. Silicon nitride ( $\text{Si}_3\text{N}_4$ ) ceramics, known for their excellent mechanical properties (high hardness, strength, and toughness), self-lubrication, high-temperature resistance, chemical stability (resistance to acids, alkalis, and molten metals), as well as transparency and wave-transmitting capabilities, are widely used in mechanical, automotive, aerospace, biomedical, and electronic applications, such as cutting tools, ceramic bearings, turbine rotors, and heat-dissipating substrates.<sup>(1-3)</sup>

During preparation, the purity of raw materials, formulation systems, and sintering methods critically influence lattice arrangement, bulk oxygen content, phase transformation, self-diffusion, and densification, thereby determining ceramic quality.<sup>(4)</sup> Producing high-performance  $\text{Si}_3\text{N}_4$  ceramics must start with high-quality  $\text{Si}_3\text{N}_4$  powder.

$\text{Si}_3\text{N}_4$  exists in three crystalline phases:  $\alpha$ ,  $\beta$ , and  $\gamma$ . The  $\alpha$ -phase, with low symmetry, high internal strain, and elevated free energy, is thermodynamically metastable.<sup>(5)</sup> Irreversible  $\alpha \rightarrow \beta$  phase transformation occurs at high temperatures (1400–1800°C). The  $\gamma$ -phase can only be synthesized under high pressure and temperature. Thus,  $\alpha$ -phase  $\text{Si}_3\text{N}_4$  dominates as the primary raw material.

Current synthesis methods for  $\alpha$ -phase  $\text{Si}_3\text{N}_4$  powder include the silicon powder nitridation method<sup>(6)</sup> and ammonolysis.<sup>(7)</sup> Both methods can produce highly pure  $\alpha$ -phase  $\text{Si}_3\text{N}_4$  materials, and the physical properties can be tuned by various synthesis parameters. For high-performance ceramics raw powders must exhibit large specific surface area and small particle size. A larger specific surface area enhances self-diffusion during sintering, promoting  $\alpha$ -phase particle rearrangement, dissolution-precipitation ( $\alpha \rightarrow \beta$  phase transformation), and solid-state diffusion, leading to uniform microstructures. Therefore, it is essential to characterize the specific surface area of powders synthesized via methods such as silicon powder nitridation and ammonolysis using the  $\text{N}_2$  physical adsorption method.

## **2. Experiment**

BET specific surface areas were calculated with  $\text{N}_2$  adsorption-desorption experiments conducted on commercial silicon nitride ( $\text{Si}_3\text{N}_4$ ) powders with varying specific surface areas, denoted A-H. The first round of tests (A-F) used a static physisorption analyzer from **AMI**, and the second round of tests (G-H) used the **AMI Surface DX 400**, a dynamic flow adsorption analyzer invented by **AMI**.

### 3. Results

The results from the static N<sub>2</sub> physisorption analyzer are summarized in Table 1. The relative standard deviation (RSD) for all six samples was less than 1%, which fully meets industrial requirements of internal standards for specific surface area measurements.

Si <sub>3</sub> N <sub>4</sub> Sample	Specific Surface Area (m <sup>2</sup> /g)			Average (m <sup>2</sup> /g)	Standard Deviation (m <sup>2</sup> /g)	RSD (%)
	Test Number					
	1	2	3			
A	2.609	2.587	2.594	2.597	0.009	0.35
B	4.522	4.578	4.569	4.556	0.025	0.54
C	4.596	4.676	4.661	4.644	0.035	0.75
D	2.620	2.660	2.669	2.650	0.021	0.80
E	2.277	2.299	2.298	2.291	0.01	0.44
F	0.429	0.435	0.427	0.430	0.003	0.79

Table 1: Experimentally determined specific surface areas for Si<sub>3</sub>N<sub>4</sub> samples A-F using a traditional static physisorption analyzer

For comparison, commercial Si<sub>3</sub>N<sub>4</sub> samples G and H were tested with the **AMI Surface DX 400** dynamic flow system, and the results are shown in Table 2.

Si <sub>3</sub> N <sub>4</sub> Sample	Specific Surface Area (m <sup>2</sup> /g)						Average (m <sup>2</sup> /g)	Standard Deviation (m <sup>2</sup> /g)	RSD (%)
	Test Number								
	1	2	3	4	5	6			
G	9.164	9.222	9.213	9.158	9.162	9.218	9.190	0.031	0.34
H	1.840	1.846	1.852	1.847	1.848	1.839	1.845	0.005	0.27

Table 2: Experimentally determined specific surface areas for Si<sub>3</sub>N<sub>4</sub> samples G-H using the dynamic flow adsorption analyzer (**AMI Surface DX 400**)

The calculated specific surface areas from the AMI Surface DX had very low standard deviation and RSD values, matching the precision of the traditional physisorption analyzer. However, the AMI Surface DX uses a dynamic flow method that can calculate single point and multipoint BET specific surface area. The unique dynamic flow method measures adsorption much more rapidly, up to 8 samples an hour, while static BET analysis often requires several hours for a single isotherm.

## 4. Conclusions

To enhance testing efficiency in industrial production lines, **AMI** has developed the high-efficiency dynamic volumetric method (**Surface DX series**) for rapid detection of specific surface area. Replicate tests on silicon nitride ( $\text{Si}_3\text{N}_4$ ) powders with varying specific surface areas demonstrate a relative standard deviation (RSD) of less than 1%, fully meeting the requirements for efficient quality control of in-line products for industrial clients. This instrument is highlighted in Figure 1.




Figure 1: Highlighted features of the **AMI Surface DX 400**

## 5. References

- (1) Du, X.; Lee, S. S.; Blugan, G.; Ferguson, S. J. Silicon nitride as a biomedical material: An overview. *Int. J. Mol. Sci.* **2022**, *23*, 6551.
- (2) Nakashima, Y.; Miyazaki, H.; Zhou, Y.; Hirao, K.; Ohji, T.; Fukushima, M. Sintered reaction-bonded silicon nitride ceramics for power-device substrates - Review -. *Open Ceram.* **2023**, *16*, 100506.
- (3) Kandi, K. K.; Thallapalli, N.; Chilakalapalli, S. P. R. Development of silicon nitride-based ceramic radomes – A review. *Int. J. Appl. Ceram. Technol.* **2014**, *12*, 909-920.
- (4) Heimann, R. B. Silicon nitride ceramics: Structure, synthesis, properties, and biomedical applications. *Materials*, **2023**, *16*, 5142.
- (5) Goto, Y. and Thomas, G. Phase-transformation and microstructural changes of  $\text{Si}_3\text{N}_4$  during sintering. *J. Mater. Sci.*, **1995**, *30*, 2194-2200.
- (6) Chang, F.-W.; Liou, T.-H.; Tsai, F.-M. The nitridation kinetics of silicon powder compacts. *Thermochim. Acta*, **2000**, *354*, 71-80.
- (7) Mazdiyasi, K. S. and Cooke, C. M. Synthesis, characterization, and consolidation of  $\text{Si}_3\text{N}_4$  obtained from ammonolysis of  $\text{SiCl}_4$ . *J. Am. Ceram. Soc.* **1973**, *56*, 628-633.